


# Use of Innovative Information Systems Combining HIV-1 Genotypic and Phenotypic Drug Resistance Interpretations For Routine and Research Applications



**30<sup>th</sup> Annual Clinical Virology Symposium**  
 April 27 – 30, 2014 | Daytona Beach, Florida

Isabelle Hardy<sup>1</sup>, Daniela Moisi<sup>2</sup>, Hugues Charest<sup>3</sup>, Régis Cantin<sup>3</sup>, Michel Roger<sup>1</sup>, Mark A. Wainberg<sup>2</sup>, Cécile Tremblay<sup>1,3</sup>, Matthieu Barralon<sup>4</sup>, Ronan Boulmé<sup>4</sup>, Chalom Sayada<sup>4</sup>, Dimitri Gonzalez<sup>5</sup>, Bluma Brenner<sup>2</sup>

<sup>1</sup>Centre Hospitalier de l'Université de Montréal, Montreal, Quebec

<sup>2</sup>Jewish General Hospital, Montreal, Quebec


<sup>3</sup>Laboratoire de santé publique du Québec/INSPQ, Montreal, Quebec

<sup>4</sup>ABL SA, Luxembourg, Luxembourg

<sup>5</sup>ABL TherapyEdge Spain SL, Barcelona, Spain


## Background

- Long-term management of HIV-1 infection requires individualized strategies tailored to each patient profile.
- The selection of antiretroviral (ARV) regimens needs to take into account the drug resistance patterns harbored by the individual.
- We have evaluated the ViroScore-HIV drug resistance algorithmic platform for routine clinical and research activities.
- This platform integrates genotypic analysis with multiple interpretative genotypic and phenotypic algorithms (Fig. 1A) to effectively track emergent resistance profiles that confer resistance to the 6 possible drug classes (Fig. 1B).



## Methods


- ViroScore-HIV is a standardized and secured web solution, certified for In-Vitro Diagnostics (IVD) use which was installed and evaluated for the management of the routine clinical and research activities of 3 genotyping laboratories in Quebec, Canada (Fig. 2).
- Preliminary analyses on 50 sequences from our clinical database evaluated ViroScore-HIV on samples previously tested by the discontinued VircoType platform.
- The ViroScore-HIV platform was very flexible, providing algorithms for up to 8 different drug resistance guidelines spanning every protein of interest (protease, reverse transcriptase, integrase, GP41 and GP120 including tropism determination).
- Results (subtyping, variants calling and drug resistance) were compared to VircoType to assess the reliability and the user-friendliness of the application for diagnostics use.
- In addition, software was combined with VisibleChek, a data mining application tailored for the management of Sanger and Next Generation Sequencing (NGS) data, to assess usability for future research and clinical applications (Fig. 4).



**Fig. 2:** Overview of the ViroScore installation and IT infrastructure.

## Results

- ViroScore-HIV was installed on a central server and configured according to the needs of each site to optimize data workflows and improve timing constraints.
- During November 2013, 50 samples coming from our reference laboratories were first analyzed.
- In most cases, Virtual-phenotyping interpretations from the latest version of the geno2pheno algorithm recently integrated in ViroScore gave results comparable to the VircoType interpretations.
- The provincial ViroScore-HIV genotyping reports provides a virtual-phenotyping algorithm (geno2pheno[resistance] v3.3) with mutational interpretations from the Stanford HIVdb v6.3.1 algorithm.
- The ViroScore-HIV subtyping tool has been updated to include a recently developed High Performance Subtyping tool.
- As the intent is to use ViroScore-HIV for daily diagnostic use, reports (Fig. 3) were exchanged with clinicians who confirmed their usefulness.
- Since January 2014, ViroScore-HIV has been used in routine by our laboratories to generate 115 reports.



**Fig. 3:** Diagnostic use configuration of the ViroScore-HIV report.



**Fig. 4:** Data mining through VisibleChek for research and clinical studies applications.

## Conclusions

- Using the advanced ViroScore-HIV information systems, combining genotypic and phenotypic algorithms for drug resistance determinations helped in routine diagnostics and facilitated the daily management of HIV infection.
- In addition, analysis of the results within VisibleChek, may assist in defining emergent resistance patterns for newer second and third generation antiretroviral drugs.

- Furthermore, the technology also enables research activities via the use of easy-to-use interfaces and data workflows tailored to the specificity of a network of hospitals.