

Evaluation of a commercial assay whole genome HIV-1 using next-generation sequencing



for the detection of HIV-1 drug resistance mutations

O. Ardizzoni¹, L. Naldi¹, J. Bernard¹, L. Deblir², D. Gonzalez¹, R. Boulme¹, C. Sayada², S. Mohamed¹.

¹ABL FRANCE - Marseille (France), ²ABL LUXEMBOURG - Luxembourg (Luxembourg)


32nd **ECCMID** Lisbon, Portugal
23–26 April 2022

Background

Drug-resistance mutations are routinely detected using standard Sanger sequencing, which does not detect minor variants with a frequency below 15%. NGS is thus becoming the new standard for genotypic drug resistance testing for HIV. The global spread of SARS-CoV-2 mobilized both the public and private sector and resulted in a rapid development of solutions focused on SARS-CoV-2 detection and sequencing. Many laboratories are now equipped to perform Whole Genome Sequencing (WGS). The objective of this study was to evaluate the performance of the WGS of HIV-1 assay.

Methods

A total of 10 positive samples were prepared, extracted (Roche MagNa Pure 24, Roche). Three HIV-1 genomic targets were amplified using the CEIVD DeepChek® Assays PR/RTand INT regions (ABL) and were sequenced using the SANGER SeqStudio system (Applied Biosystem). The whole genome HIV-1 was amplified using the DeepChek® Assays Whole Genome HIV-1 (ABL) and was sequenced using the NGS iSeq100 (Illumina). Sequences were compared to those obtained by Sanger Sequencing. HIV-1 QCMD and AcroMetrix™ HIV Mutant (Thermofisher scientific) were sequenced. DeepChek® HIV-1 Whole Genome software (ABL) was used for the interpretation of drug resistance.



Results

The median coverage per sample for the WGS of HIV-1 was 17'500 reads. High analytical reproducibility and repeatability were evidenced by Percent Agreement being 100%. Duplicated samples in two different NGS runs were 100% homologous. NGS detected all the mutations found by Sanger sequencing and identified additional resistance variants. The score of the QCMD panel detection of drug resistance mutations for RT/PR and INT were 339/339 and 125/125, respectively. All AcroMetrix™, Sersesq mutations were detected using WG HIV kit

Sample	Ref	Viral Load	subtype	RT mutations of interest	PR mutations of interest	INT mutations of interest
1	HIV Zeptomatrix	2.10 ⁷ TCID50	B	wildtype	wildtype	wildtype
2	HIV Zeptomatrix	2.10 ⁵ TCID50	B	wildtype	wildtype	wildtype
3	HIV Zeptomatrix	2.10 ¹ TCID50	B	wildtype	wildtype	wildtype
4	HIV Acromatrix	140,000 cp/mL	B	M41L, T69D, K70R, Y181C, T215Y	L10R, D30N, A71V, L90M	No mutation
5	HIV Acromatrix	28,000 cp/mL	B	M41L, T69D, K70R, Y181C, T215Y	L10R, D30N, A71V, L90M	No mutation
6	HIV Acromatrix	20,000 cp/mL	B	M41L, T69D, K70R, Y181C, T215Y	L10R, D30N, A71V, L90M	No mutation
7	HIV Sersesq	20,000cp/mL	B	L10I	T215Y, K219Q	No mutation
8	QCMD1	4.57 Log10 cp/mL	C	M41L, E44D, D67N, T69D, A98G, M184V, L210W, T215Y	L10F, D30N, N88D	No mutation
9	QCMD2	4.24 Log10 cp/mL	A/G	No mutation	No mutation	No mutation
10	QCMD3	5.16 Log10 cp/mL	B	No mutation	K43T (19,63%), M46I, I54V, V82A, L90M	No mutation

Conclusions

This study is the first evaluation of the DeepChek® Assays Whole Genome HIV-1 (ABL) using the iSeq100 system combined with an easy software. The NGS should occupy a major place in HIV resistance surveillance and clinical care, thanks to its decreasing costs (due to COVID-19 pandemic) and ability to reveal resistant variants and study their impact; especially on the new capsid/maturation inhibitors and detection of potential new clinically relevant mutations in the HIV genome.

References

- Mohamed S, Penaranda G, Gonzalez D, Camus C, Khiri H, Boulme R, et al. Comparison of ultra-deep versus Sanger sequencing detection of minority mutations on the HIV-1 drug resistance interpretations after virological failure. AIDS 2014,28:1315-1324.
- Mohamed S, Ravet S, Camus C, Khiri H, Olive D, Halfon P. Clinical and analytical relevance of NNRTIs minority mutations on viral failure in HIV-1 infected patients. J Med Virol 2013; 86:394–403.

Contact information

Sofiane Mohamed, PhD:
s.mohamed@ablsa.com

<https://www.ablsa.com/>