




RapidChek SARS-CoV-2 Antibody Rapid Test Kit



“COLLOIDAL GOLD IMMUNOCHROMATOGRAPHY”



ACCURATE

- High purity antibody with high accuracy



EFFICIENT

- Result available within 15 minutes
- Shelf life 12 months



CONVENIENT

- Suitable for finger blood.
- Storage at 4°C – 30°C



Imported by **Advanced Biological Laboratories (ABL) S.A.**

17 rue des Jardiniers, L-1835 Luxembourg

Technical Support Center at support-diag@ablsa.com



Call +339 7017 0300 or contact your ABL Field-Application Specialist




Kit content

Each kit contains 20 tests, for each test :

- One testing strip
 - Gold Standard mat (colloidal gold labeled COVID-19 recombinant protein)
 - Sample mat
 - Cellulose nitrate membrane IgM & IgG
- One dropper
- One package of desiccant



- C: Quality control
- G : IgG detection line
- M : IgM detection line



References & contact



RapidChek SARS-CoV-2 Antibody Rapid Test Kit. K-20-RC-CoV-2

CONTACT

- Advanced Biological Laboratories S.A.
- 17, rue des Jardiniers, L-1835 Luxembourg
- contact@ablsa.com
- <http://www.ablsa.com>
- TEL : (+352) 26 38 96 761
- FAX : (+352) 26 38 96 76 31



Imported by **Advanced Biological Laboratories (ABL) S.A.**

17 rue des Jardiniers, L-1835 Luxembourg

Technical Support Center at support-diag@ablsa.com

Call +339 7017 0300 or contact your ABL Field-Application Specialist